

EEG PREPROCESSING

1. what is Preprocessing

Preprocessing refers to the set of transformations applied between raw EEG data collection and actual analysis. The goal is to maximize the **signal-to-noise ratio** while preserving meaningful brain signals.

2. Preprocessing Pipeline

Based on the textbook chapter, the recommended order of EEG preprocessing is:

1. Filtering (Continuous Data Only)

Removes noise such as slow drifts (high-pass filter $\sim 0.1\text{--}0.5$ Hz), muscle artifacts (low-pass filter), and electrical line noise (50/60 Hz notch filter).

- **High-pass filter ($\sim 0.1\text{--}0.5$ Hz)** → removes slow drifts.
- **Low-pass filter (40–80 Hz)** → removes muscle and high-frequency noise.
- **Notch filter (50/60 Hz)** → removes electrical line noise.

Why First?

Filtering on continuous data avoids edge artifacts that occur if filtering is applied after epoching.

2. Re-referencing

- EEG is measured relative to a reference electrode.
- Common choices: linked mastoids, earlobes, or average reference.
- This reduces global artifacts and makes signals more spatially specific.
- Poor referencing can bias signals across the scalp.

Why Here ?

- Referencing is done early so that subsequent steps (epoching, artifact rejection) are applied to clean, referenced data.

3. Epoching

- Continuous EEG is segmented into time-locked trials (stimulus onset or response).
- Epoch length must cover both pre- and post-event activity.
- Include a **baseline window** (e.g., -200 to 0 ms).
- **Epoch length** depends on analysis:
 - ERP → short epochs (-200 to 800 ms).
 - Time-frequency → longer epochs with buffer zones to avoid edge effects.

Why After Referencing?

- Trials are cleaner and more consistent when referenced before cutting into epochs.

4. Baseline Correction

- Subtract the mean activity in the pre-stimulus baseline window from the entire epoch.
- **Purpose:**
 - Removes slow drifts.
 - Normalizes trials for fair comparison.

Why After Epoching?

- Baseline correction requires a per-trial reference (pre-stimulus period).

5. Artifact Rejection / Correction

- Methods:
 - Trial rejection (automatic/manual removal of noisy trials).
 - ICA → removes eye blinks, muscle artifacts while preserving brain signals.

Why Here?

- Perform artifact rejection **after epoching and baseline correction**, so that trials/channels with strong artifacts do not contaminate averages.

6. Interpolating Bad Electrodes

- Replace noisy/dead channels with interpolated values from neighbors.
- Better to fix during recording

Why Here?

Ensures no channel artificially biases averages.

7. Spatial Filtering(Optional)

- Examples: Surface Laplacian..
- **Purpose:** Improve localization and reduce volume conduction.

Why Hater?

- Applied only once data are relatively clean, often just before time-frequency or connectivity analysis.

8. Trial Count Matching (Optional, for experiments with conditions)

- Important for condition comparisons: ensures one condition doesn't appear stronger simply because it had more trials.

- Equalize number of trials across conditions (random selection, first-N, or based on behavioural metrics like reaction time etc).

Why Here?

- Prevents bias: conditions with more trials shouldn't appear stronger just due to averaging.